



ICMGP 2024
CAPE TOWN • SOUTH AFRICA • 21 - 26 JULY

Seth Lyman, Tyler Elgiar, Utah State University
**Advancing Permeation Tube-based Mercury
Calibration**



Acknowledgements

- Funding provided by the U.S. National Science Foundation
- This work involves and owes much to:
 - Utah State University: Trevor O'Neil, KarLee Zager
 - Reed College: Lynne Gratz, Ellie Derry
 - University of Nevada, Reno: Mae Gustin, Sarrah Dunham-Cheatham, Jan Gacnik (now at Jožef Stefan Institute), Livia Lown
 - Jožef Stefan Institute: Milena Horvat, Igor Zivkovic, Sreekanth Vijayakumaran Nair, Teodor Andron
 - University of Colorado Boulder: Rainer Volkamer
 - University of Utah: Gannet Hallar

(Tyler Elgiar is now with the U.S. Bureau of Land Management)



The Problem



The Problem

- A. Commercial (i.e., KCl denuder-based) measurement systems for Hg^{II} are, at best, merely qualitative



The Problem

- A. Commercial (i.e., KCl denuder-based) measurement systems for Hg^{II} are, at best, merely qualitative
- B. It took a decade and more to figure this out because no field-deployable calibration system existed



The Problem


- A. Commercial (i.e., KCl denuder-based) measurement systems for Hg^{II} are, at best, merely qualitative
- B. It took a decade and more to figure this out because no field-deployable calibration system existed
- C. We need methods for SI-traceable Hg^{II} calibration so we don't repeat past mistakes



Options for SI-traceable calibration

Options for SI-traceable calibration

1. Aspiration or evaporation of liquid standards



Sari et al., 2019, Gacnik et al., 2021



Options for SI-traceable calibration

2. Flow over a saturated bed




PS Analytical Mercury Vapour Generator



Options for SI-traceable calibration

3. Diffusion tubes




Xiao et al., 1997; Swartzendruber et al., 2009; de Krom et al. (multiple; for Hg⁰)



Options for SI-traceable calibration

4. Direct Hg^0 oxidation




Gacnik et al., 2022



Options for SI-traceable calibration

5. Permeation tubes



First used by Landis et al., 2002?

Why are we putting our effort into permeation tube methods?



Why are we putting our effort into permeation tube methods?


- Capable of low, continuous emission rates
- Shown to work in real field conditions
- Confidence in identification of compound emitted
- Simple (in principle)



Basic principles of our method




Basic principles of our method




Basic principles of our method

- Tubes kept warm (50-70°C?), which keeps the permeation rate stable
- Permeable area minimized to keep permeation rate low
- Inert materials
- Give ample time for equilibration





Verification of chemical speciation by GCMS





Recent success (Elgiar et al., 2024): Validation of SI-traceability



Recent success (Elgiar et al., 2024): Validation of SI-traceability




Recent success (Elgiar et al., 2024): Validation of SI-traceability



Remaining Challenge: All available materials retain Hg compounds


Remaining Challenge: All available materials retain Hg compounds



- HgBr₂ is retained by all surfaces, even at high temperature (up to 150°C)
- About 2 hours of continuous flow is required to achieve equilibration

Remaining Challenge: The last centimeter


Remaining Challenge: The last centimeter



- HgBr_2 can be retained on unheated inlet tip
- Hg^{II} can spike when inlet is first inserted

Remaining Challenge: Stability/permeation rate tradeoff


Remaining Challenge: Stability/permeation rate tradeoff



- Permeation rates are:
 - unstable at low temperatures (probably even at 50°C)
 - Too high at high temperatures
- Intermittent injection is difficult
- Dilution calibrators could work

Remaining Challenge: How to make speciation of output NIST-traceable?

Remaining Challenge: How to make speciation of output NIST-traceable?



- Our method results in some Hg⁰ output
- Is there a traceable method to verify how much?

Remaining Challenge: Will it work for non-halide Hg compounds?



