



ICMGP 2024
CAPE TOWN • SOUTH AFRICA • 21 - 26 JULY
CAPE TOWN INTERNATIONAL CONVENTION CENTRE

Stability of Common Preconcentration Methods for Gaseous Oxidized Mercury in Air

Sreekanth Vijayakumaran Nair, Saeed Waqar Ali, Jan Gačnik, Igor Živković,
Teodor Daniel Andron, and Milena Horvat*

Jožef Stefan Institute, Jamova cesta 39, Ljubljana - Slovenia



INTRODUCTION

Atmospheric mercury (Hg): Elemental form (**GEM**)

Oxidized forms (Particle bound, **PBM** & Gaseous oxidized, **GOM**)

Ultra-trace level concentration and high reactivity

GOM measurements: requires a pre-concentration step

Biased low and lack of comparable measurements results

- Biases from improper calibration (traceability)
- Biases from sampling methods

Accurate measurements – critical for Minamata Convention effectiveness evaluation!!!


Common methods

Denuders (1 – 2 hours)

&

CEMs

(days to weeks)



under-estimation
of Hg^{II} in air

- used for sampling GOM (Hg^{II})
- sampling losses → **biases**

Denuders: analysis by thermal decomposition of Hg^{II} to Hg^0 and detection by CVAFS

CEMs: Overnight digestion in BrCl solution, followed by analysis using CVAFS

OBJECTIVES

➤ Laboratory experiments were performed using ^{197}Hg radiotracer




Advantages:

- Enhanced experimental precision
- No blank or contamination issues
- High specificity and ease of detection
- Experiments at environmentally relevant concentration


➤ Quantification of biases: based on precise mass balance budget

- Evaluate the performance of denuders & CEMs in retaining different GOM species
- Assess the recovery efficiency of GOM from denuders and CEMs


METHODOLOGY



A) Losses from denuders



B) Losses from CEM




Denuders

1. Hg^{II} losses captured on Au & KCl traps (during exposure)
2. Hg^{II} recovered from denuders (post-exposure)




CEMs

1. Hg^{II} losses captured on Au & KCl traps (during exposure)
2. Hg^{II} recovered by digestion (post-exposure)
3. Residual Hg^{II} recovered from CEM (post-digestion)
4. Hg^{II} recovered from inner Teflon parts (acid-wash)

Plasma loading efficiency - denuders






- NTP loading efficiency ranged between 86.2–96.1%.
Avg. efficiency:
DEN A – 94.2%
DEN B – 89.0%
DEN C – 90.0%
- Avg. breakthrough:
DEN A – 6.5%
DEN B – 10.6%
DEN C – 10.4%
- GOM recovery by thermal decomposition (500 °C ~ 15 min)
DEN A – 101%
DEN B – 100%
DEN C – 99.9%



➤ Avg. losses 30.1%, 35.9%, and 35.1% for denuder A, B, and C respectively

Condition: Ambient air (dark)




➤ Avg. losses 19.5%, 19.1%, and 25.5% for denuder A, B, and C respectively

Condition: Ambient air (light exposure)

Summary (denuders)

- Freshly prepared denuders lost upto **36% GOM (Hg^{II} reduced to Hg^0)** <<< when exposed to ambient air (exposure period – 2 hours)
 - No losses were observed in the form of Hg^{II} (no Hg^{II} captured on KCl traps)
 - More humidity in air, more losses
 - Minimum **64%** of Hg^{II} loaded in denuders were left for measurement after the exposure period
- Re-used denuders lost upto **79.7% GOM**; Avg. loss was 51%

Plasma loading efficiency - CEMs



NTP oxidation efficiency

- 90.2–102% for HgO
- 94.2–98.0% for HgCl₂
- 92.2–98.2% for HgBr₂

CEM digestion (Avg. recoveries)

- 94.8% for HgO
- 94.0% for HgCl₂
- 93.7% for HgBr₂

Teflon walls (Avg. recoveries)

- 6.5% for HgO
- 4.3% for HgCl₂
- 7.7% for HgBr₂

Postdigestion residuals on CEM (Avg.)

- 5.5% for HgO
- 6.4% for HgCl₂
- 6.7% for HgBr₂

i) Hg^{II} losses from CEM (72 hours): HgO loaded on CEM

Amount (ng)	Losses (%)	CEM Digestate (%)	CEM Leftover (%)	FP wash (%)	Mass balance (%)
2.5/60	17.6	64.8	7.5	5.0	94.8
10	18.3	74.4	5.7	1.4	99.9
30	14.4	71.8	6.4	6.2	98.8
Average	16.7	70.3	6.5	4.2	97.8

ii) Hg^{II} losses from CEM (72 hours): HgCl₂ loaded on CEM


Amount (ng)	Losses (%)	CEM Digestate (%)	CEM Leftover (%)	FP wash (%)	Mass balance (%)
2.5	3.1	86.1	5.6	3.0	97.7
10.0	3.0	87.2	4.0	2.5	96.7
30.0	2.1	90.6	5.0	1.2	99.0
Average	2.7	88.0	4.9	2.2	97.8

iii) Hg^{II} losses from CEM (72 hours): HgBr₂ loaded on CEM

Amount (ng)	Losses (%)	CEM Digestate (%)	CEM Leftover (%)	FP wash (%)	Mass balance (%)
2.5	3.1	92.4	4.7	1.0	101.2
10.0	5.2	86.2	4.0	3.6	99.0
30.0	3.1	92.4	4.7	1.0	101.2
Average	3.8	90.4	4.5	1.9	100.5

Summary (CEMs)

- CEM-based measurements demonstrate greater reliability over denuders (lower losses of Hg^{II} compounds)
 - HgO – 16.7% (average loss)
 - HgCl₂ and HgBr₂ – 2.7% and 3.8%, respectively (average loss)
- CEMs may not uniformly collect all types of GOM species with equal efficiencies
 - better performance in retaining cationic species of GOM (chlorides and bromides)
 - not optimal for neutral species such as HgO
- Hg^{II} recovery from CEMs may not be adequate (few % lost to Teflon walls and CEMs as residuals; post-digestion)
 - Alternative methods:
 - 1) Leaching using tetrabutylammonium (commonly used for aqueous biphasic extraction of metal ions) – recovery improved by 3.5%, but not complete.
 - 2) ultra-wave digestion method recovered upto 100.4% residuals from CEMs
 - future studies to assess its robustness and potential as a viable alternative to BrCl digestion



Thank you

This research work was part of the GMOS-Train project funded by European Union's Horizon 2020 research and innovation programme under the Marie Skłodowska-Curie grant agreement no. 860497. The research has also been co-funded by the SI-Hg project under the European Metrology Programme for Innovation and Research (EMPIR; grant agreement no. 19NRM03), Slovenian Research and Innovation Agency (ARIS) programme P1-0143, and projects J1-1716 (STRAP) and J1-3033 (IsoCont).